



EASY MOVE

POWER WHEELCHAIRS

HIGH QUALITY
ATTRACTIVE PRICE





EM002



Type: EM002

The basic power wheelchair



Features:

- Quickly detachable battery box
- Flip-up armrests
- Swing-away footrests
- Joystick interchangeable left and right
- Vertically foldable
- Fully electrical
- Manual pushing mode available
- Front Solid tyres
- Rear Air tyres

Specifications

Type	EM002
Carrying capacity	100 kg
Maximum speed	5-7 kmh
Length	85 cm
Width	65 cm
Motor	2x 250W
Weight with batteries	44 kg
Weight without batteries	35 kg
Maximum range	12(±5) km*
Battery	2x 12 Ah
Recharge time	6-8hrs
Maximum gradient	8-15 degrees

* Range varies by terrain, weight of rider and ambient temperature

EM23-16



Type: EM23-16

The foldable power wheelchair



Features:

- Quickly detachable battery box
- Fixed armrests
- Swing-up footrest
- Joystick interchangeable left and right
- Completely foldable
- Fully electrical
- Manual pushing mode available
- Front Solid tyres
- Rear Air tyres

Specifications

Type	EM23-16
Carrying capacity	90 kg
Maximum speed	5-7 kmh
Length	97 cm
Width	63 cm
Motor	2x 200W
Weight with batteries	29 kg
Weight without batteries	20 kg
Maximum range	12(±5) km*
Battery	2x 12 Ah
Recharge time	6-8 hrs
Maximum gradient	8-15 degrees

* Range varies by terrain, weight of rider and ambient temperature



Easy to carry, transportable power wheelchair for short and medium journeys.





EM23-41



Type: EM23-41

The foldable power wheelchair with Li-ion battery 

Features:

- Quickly detachable Lithium battery
- Flip-up armrests
- Swing-up footrest
- Joystick interchangeable left and right
- Compact foldable
- Fully electrical
- Manual pushing mode available
- Front Solid tyres
- Rear Air tyres

Specifications

Type	EM23-41
Carrying capacity	90 kg
Maximum speed	8 kmh
Length	87 cm
Width	62 cm
Motor	2x 200W
Weight with batteries	22 kg
Weight without batteries	19 kg
Maximum range	12(±5) km*
Battery	LITHIUM
Recharge time	2-3 hrs
Maximum gradient	8-15 degrees

* Range varies by terrain, weight of rider and ambient temperature

EM23-20



Type: EM23-20

The luxury power wheelchair 

Features:

- Detachable battery box
- Easy to take out armrests
- Swing-away footrests
- Joystick interchangeable left and right
- Vertically foldable
- Fully electrical
- Manual pushing mode available
- Front Air tyres
- Rear Air tyres

Specifications

Type	EM23-20
Carrying capacity	100 kg
Maximum speed	8-10 kmh
Length	97 cm
Width	65 cm
Motor	2x 320W
Weight with batteries	63 kg
Weight without batteries	40 kg
Maximum range	30(±5) km*
Battery	2x 35 Ah
Recharge time	6-8 hrs
Maximum gradient	8-15 degrees

* Range varies by terrain, weight of rider and ambient temperature



The comfort power wheelchair for both indoor and outdoor.



EM23-28L



Type: EM23-28L

The easy folding power wheelchair with Li-ion battery 

Features:

- Quickly detachable Lithium battery
- Flip-up armrests
- Swing-up footrest
- Joystick interchangeable left and right
- Compact foldable
- Fully electrical
- Manual pushing mode available
- Front Solid tyres
- Rear Solid tyres

Specifications

Type	EM23-28L
Carrying capacity	110 kg
Maximum speed	8 kmh
Length	87 cm
Width	56 cm
Motor	2x 200W
Weight with batteries	22 kg
Weight without batteries	19 kg
Maximum range	16 km*
Battery	LITHIUM
Recharge time	2-3 hrs
Maximum gradient	8-15 degrees


* Range varies by terrain, weight of rider and ambient temperature



EM30TT



Type: EM30TT

The multifunctional power wheelchair 

Features:

- Standing possibility
- Reclining possibility
- Flip-up armrests
- Detachable head rest
- Joystick interchangeable left and right
- Fully electrical
- Manual pushing mode available
- Front Solid tyres
- Rear Air tyres

Specifications

Type	EM30TT
Carrying capacity	120 kg
Maximum speed	8-10 kmh
Length	105 cm
Width	66 cm
Motor	2x 320W
Weight with batteries	82 kg
Weight without batteries	60 kg
Maximum range	30(±5) km*
Battery	2x 35 Ah
Recharge time	6-8 hrs
Maximum gradient	8-15 degrees

* Range varies by terrain, weight of rider and ambient temperature



Standing and reclining power wheelchair for both indoor and outdoor





Can we help?

We appreciate there are a lot of things to think about when selecting the right Power Wheelchair and because of that we are more than happy to answer any questions you may have. You can call our advice line or have one of our friendly trained advisors come out and visit you. To ask any questions or arrange a no obligation demonstration on location contact us today.

We are happy to help.

Easy Move: Whitefield, Bangalore, Karnataka, India

TELEPHONE +91-7406062888 or +91-7406135888
info@easymove-mobility.com | www.easymove-mobility.com



All Easy Move models are CE certified and classified as 'Electric motor driven wheelchairs with manual steering'.